

AMPLIFY[®]
SURGICAL

SEEING IS BELIEVING



dualPortal[®]



dualX[®]

A novel **two-portal endoscopic** approach to the spine

The Problem




Common Spine Surgeries Today

Less than 25% of all spine procedures are performed in a minimally invasive way, 4% Endoscopically

	COMMON SPINAL SURGERY	COMMON MIS SURGERY
BLOOD LOSS	1018 mL	200mL
INTRA-OP COMPLICATIONS	9%	1.6%
HOSPITAL STAY	5 days	< 24 hours
INCISION LENGTH	10cm	< 1cm





Endoscopy in Other Surgeries



Endoscopy in Spine - Challenges

- Too difficult
- Too expensive
- Limited to only a small % of procedural types
- No robust solutions for more common procedures
 - e.g. Lumbar Fusions
- Inferior surgical implant options for endoscopic fusion



The Solution

dualPortal Endoscopic Spine

dualPortal[®]




THE dualPortal[®] SPINAL ENDOSCOPY

The dualPortal[®] solution is a novel two-portal endoscopic approach to the spine that allows surgeons to easily learn and perform a wider array of lumbar spine procedures with lower cost compared to the conventional one-portal technique. It also provides flexibility to perform endoscopic lumbar fusions with the dualX[®] Expanding Interbody Fusion System.


Why dualPortal[®] Spine?

Compared to uni-portal spine endoscopy.

- **Shorter Learning Curve:** similar workflow to a familiar posterior approach.
- **Lower Cost:** compatible with widely available endoscopic systems.
- **Versatility:** accommodates a variety of lumbar procedures, including multi-level fusions.



dualPortal[®]




+ dualX[®]



Greater
Flexibility



Enhanced
Visualization



Shallow
Learning Curve

dualPortal Access System

AMPLIFY
SURGICAL

dualPortal®



THE dualPortal® SPINAL ENDOSCOPY

The dualPortal® solution is a novel two-portal endoscopic approach to the spine that allows surgeons to easily learn and perform a wider array of lumbar spine procedures with lower cost compared to the conventional one-portal technique. It also provides flexibility to perform endoscopic lumbar fusions with the dualX® Expanding Interbody Fusion System.



dualPortal® CANNULA

- Manages irrigation fluid control.
- Critical for visualization.
- Facilitates easy access for working instruments.




dualPortal® HEMOSTATIC TUBE

Delivers hemostatic agent to targeted site.




dualPortal® SCOPE RETRACTOR


Transforms the endoscope additionally into a retractor.



TRAY 1



TRAY 2




TRAY 3

Introducing
dualPortal[®] 2.0
ENDOSCOPIC SYSTEM

DESIGNED TO BE SEEN


A COMPREHENSIVE SOLUTION FROM O.R. PREP TO FINAL CLOSE



AMPLIFY[®]
SURGICAL

DEPTH GAUGE

Measure the tissue depth
from skin to lamina.



dualPortal® 2.0

MULTI-PURPOSE GUIDE

Provide controlled guidance on proper position, angle, and depth for creating the two incisions needed for dualPortal® surgery

Provide guidance and safety stop


Provide guidance for the scope and instruments to triangulate and "meet" at the correct site.

4~5mm incision through skin for the Scope

9~10mm incision through muscle fascia for the Working Portal

Target Needle

dualPortal® 2.0




ANGLE GUIDE

Provide controlled guidance on proper position, angle, and depth for creating the two incisions needed for dualPortal® surgery



TISSUE DISSECTOR

Dissect and loosen soft tissue to prepare space for creation of the endoscopic working space





ERGONOMIC SCOPE HANDLE

Ergonomic Scope Handle to
distribute weight and provide
support


Detachable Stabilizer Legs for
additional support


dualPortal® 2.0

ERGONOMIC SCOPE HANDLE

Ergonomic Scope Handle to
distribute weight and provide
support

Detachable Stabilizer Legs
for additional support





dualPortal® 2.0 CANNULA

Cannula provides retraction at the incision site for easier instrument insertion.

dualPortal® 2.0

DUALPORTAL 2.0 CANNULA


Cannula provides retraction at the incision site for easier instrument insertion

Active in-flow and out flow for even better visualization

Expanding tip for secure placement & retention.

Distal flare for additional soft tissue retraction

dualPortal® 2.0



A 3D digital rendering of a surgical navigation system. The background is a detailed, light-colored model of a human skull. In the center, a blue surgical navigation frame is mounted on the skull. A vertical instrument with a black handle and a white tip is positioned within the frame. To the right, a pair of black surgical forceps is shown in a suspended position. In the lower right, another blue frame component is visible. The scene is lit with a soft, blueish light, creating a clinical and futuristic atmosphere.

FLAG GUIDE

Designed to provide additional aid in guiding instruments throughout the case and increase in efficiency.


dualPortal[®] 2.0

FLAG GUIDE

Designed to provide additional aid in guiding instruments throughout the case and increase in efficiency.




dualPortal[®] 2.0



OUTFLOW DAM

Control/Contain irrigation
overflow during procedure

dualPortal[®] 2.0



OUTFLOW DAM

Control/Contain irrigation
overflow during procedure

BONE WAX DISPENSER

A syringe-type instrument to dispense small amount of bone wax at the tip and apply it to bleeding bone.

Unique bone conforming tip.



dualPortal 2.0

MULTI-PURPOSE GUIDE

Provide controlled guidance on proper position, angle, and depth for creating the two incisions needed for dualPortal® surgery

Provide guidance and safety stop

Provide guidance for the scope and instruments to triangulate and "meet" at the correct site.




DUALPORTAL 2.0 CANNULA

Cannula provides retraction at the incision site for easier instrument insertion

Active in-flow and out flow for even better visualization

Distal flare for additional soft tissue retraction

Expanding tip for secure placement & retention.



OUTFLOW DAM

Control/Contain irrigation overflow during procedure



DEPTH GAGE

Measure the tissue depth from skin to lamina.



TISSUE DISSECTOR


Dissect and loosen soft tissue to prepare space for creation of the endoscopic working space



ERGONOMIC SCOPE HANDLE

Ergonomic Scope Handle to distribute weight and provide support

Detachable Stabilizer Legs for additional support



FLAG GUIDE

Designed to provide additional aid in guiding instruments throughout the case and increase in efficiency.



BONE WAX DISPENSER

A syringe-type instrument to dispense small amount of bone wax at the tip and apply it to bleeding bone.

Unique bone conforming tip.

Patented Procedural Solution

(19) United States (12) Patent Application Publication Choi et al.


(10) Pub. No.: **US 2024/0390031 A1**
(43) Pub. Date: **Nov. 28, 2024**

(54) **MULTI-PORTAL SURGICAL TOOLS AND SYSTEMS**
(52) U.S. CL. **A61B 173423 (2013.01); A61B 9850 CPC** (2016.02); **A61B 20173443 (2013.01); A61B 20173447 (2013.01)**
(71) Applicant: **Amplify Surgical, Inc.**, Irvine, CA (US)
(72) Inventors: **Andy Wonyong Choi**, Irvine, CA (US); **Clark Hutton**, Carlsbad, CA (US)
(21) Appl. No.: **18074784**
(22) Filed: **May 24, 2024**

Related U.S. Application Data
(60) Provisional application No. 63/504,248, filed on May 25, 2023; provisional application No. 63/611,874, filed on Dec. 19, 2023.

Publication Classification
(51) **Int. Cl.** **A61B 1734** (2006.01)
A61B 9850 (2006.01)

ABSTRACT
A multi-portal method for treating a subject's spine includes distracting adjacent vertebrae using a distraction instrument positioned at a first entrance along the subject to enlarge an intervertebral space between the adjacent vertebrae. An interbody fusion implant can be delivered into the enlarged intervertebral space. The interbody fusion implant can be positioned directly between vertebral bodies of the adjacent vertebrae while endoscopically viewing the interbody fusion implant using an endoscopic instrument. The patient's spine can be visualized using endoscopic techniques to view, for example, the spine, tissue, instruments, and implants before, during, and after implantation, or the like. The visualization can help a physician throughout the surgical procedure to improve patient outcome.



(10) Patent No.: **US 11,464,648 B2** (45) Date of Patent: **Oct. 11, 2022**

(56) **References Cited**
U.S. PATENT DOCUMENTS
5,171,279 A * 12/1992 Mathews (A61B 2445) 128,988
5,390,683 A * 2/1995 Pokuandi (A61B 2445) 128,988
5,762,629 A 6/1998 Kambara (A61B 2445) 128,988
5,782,832 A * 7/1998 Larson et al. (A61B 2445) 128,988
(Continued)

FOREIGN PATENT DOCUMENTS
CN 105636555 A 6/2016 (Continued)
EP 3016617 A2 5/2016 (Continued)

OTHER PUBLICATIONS
Choi, Chang Myung et al. "How I do it? Bipodal endoscopic spinal surgery (BESS) for treatment of lumbar spinal stenosis." *Acta Neurochir* (2016) 158:459-463, published Jan. 18, 2016.
(Continued)

Prior Publication Data
US 2021/0068975 A1 Mar. 11, 2021

Int. Cl. **A61B 2444** (2006.01)
A61B 1702 (2006.01)
A61B 100 (2006.01)
A61B 1768 (2006.01)

U.S. Cl. **A61F 24455** (2013.01); **A61B 10065** (2013.01); **A61B 17025** (2013.01); **A61B 20170256** (2013.01); **A61B 2017681** (2013.01)

Field of Classification Search
CPC **A61F 2444-447; A61F 24611**
See application file for complete search history.

ABSTRACT
A multi-portal method for treating a subject's spine includes distracting adjacent vertebrae using a distraction instrument positioned at a first entrance along the subject to enlarge an intervertebral space between the adjacent vertebrae. An interbody fusion implant can be delivered into the enlarged intervertebral space. The interbody fusion implant can be positioned directly between vertebral bodies of the adjacent vertebrae while endoscopically viewing the interbody fusion implant using an endoscopic instrument. The patient's spine can be visualized using endoscopic techniques to view, for example, the spine, tissue, instruments, and implants before, during, and after implantation, or the like. The visualization can help a physician throughout the surgical procedure to improve patient outcome.
(Continued)

(12) United States Patent Choi et al.

(10) Patent No.: **US 11,678,906 B2**
(45) Date of Patent: **Jun. 20, 2023**

(54) **MULTI-PORTAL SURGICAL SYSTEMS, CANNULAS, AND RELATED TECHNOLOGIES**
(71) Applicant: **Amplify Surgical, Inc.**, Laguna Hills, CA (US)
(72) Inventors: **Andy Wonyong Choi**, Irvine, CA (US); **Dong-Hwa Heo**, Seoul (KR); **Jeffrey Roh**, Seattle, WA (US)
(73) Assignee: **Amplify Surgical, Inc.**, Irvine, CA (US)
(*) Notice: Subject to any disclaimer, the term of this patent is extended or adjusted under 35 U.S.C. 154(b) by 408 days.(21) Appl. No.: **16/687,520**
(22) Filed: **Nov. 18, 2019**

References Cited
U.S. PATENT DOCUMENTS
5,390,683 A * 2/1995 Pokuandi (A61M 25/003) 600,549
5,477,637 A * 8/1995 Lober (A61M 25/003) 600,549
(Continued)

FOREIGN PATENT DOCUMENTS
CN 105636555 A 6/2016 (Continued)
EP 3016617 A2 5/2016 (Continued)

OTHER PUBLICATIONS
Choi, Chang Myung et al. "How I do it? Bipodal endoscopic spinal surgery (BESS) for treatment of lumbar spinal stenosis." *Acta Neurochir* (2016) 158:459-463, published Jan. 18, 2016.
(Continued)

Prior Publication Data
US 2021/0068663 A1 Mar. 11, 2021

Related U.S. Application Data
(65) Continuation-in-part of application No. 16/565,403, filed on Sep. 9, 2019, now Pat. No. 11,464,648.

ABSTRACT
A multi-portal method for treating a subject's spine includes distracting adjacent vertebrae using a distraction instrument positioned at a first entrance along the subject to enlarge an intervertebral space between the adjacent vertebrae. An interbody fusion implant can be delivered into the enlarged intervertebral space. The interbody fusion implant can be positioned directly between vertebral bodies of the adjacent vertebrae while endoscopically viewing the interbody fusion implant using an endoscopic instrument. The patient's spine can be visualized using endoscopic techniques to view, for example, the spine, tissue, instruments, and implants before, during, and after implantation, or the like. The visualization can help a physician throughout the surgical procedure to improve patient outcome.
(Continued)


Primary Examiner — Eduardo C Robert
Assistant Examiner — Steven J Cotroneo
(74) *Attorney, Agent, or Firm* — Perkins Coie LLP

Int. Cl. **A61B 1734** (2006.01)
A61B 500 (2006.01)
(Continued)

U.S. Cl. **A61B 173423** (2013.01); **A61B 524** (2013.01); **A61B 50893** (2013.01);
(Continued)

Field of Classification Search
CPC **A61F 2444-447; A61F 24611; A61B 17025; A61B 20170256-0262; A61B 20173445; A61B 20173447**
See application file for complete search history.

21 Claims, 20 Drawing Sheets



(12) United States Patent Choi et al.

(10) Patent No.: **US 12,226,319 B2**
(45) Date of Patent: **Feb. 18, 2025**

(54) **MULTI-PORTAL SURGICAL SYSTEMS**
(71) Applicant: **Amplify Surgical, Inc.**, Laguna Hills, CA (US)
(72) Inventors: **Andy Wonyong Choi**, Irvine, CA (US); **Dong-Hwa Heo**, Seoul (KR); **Jeffrey Roh**, Seattle, WA (US)
(73) Assignee: **Amplify Surgical, Inc.**, Irvine, CA (US)
(*) Notice: Subject to any disclaimer, the term of this patent is extended or adjusted under 35 U.S.C. 154(b) by 0 days.(21) Appl. No.: **17/992,685**
(22) Filed: **Sep. 2, 2022**

References Cited
U.S. PATENT DOCUMENTS
3,486,505 A 12/1969 Morrison
5,171,279 A 12/1992 Mathews (Continued)

FOREIGN PATENT DOCUMENTS
CN 105636555 A 6/2016 (Continued)
EP 3016617 A2 5/2016 (Continued)

OTHER PUBLICATIONS
Choi, Chang Myung et al. "How I do it? Bipodal endoscopic spinal surgery (BESS) for treatment of lumbar spinal stenosis." *Acta Neurochir* (2016) 158:459-463, published Jan. 18, 2016.
(Continued)

Prior Publication Data
US 2023/0104335 A1 Apr. 6, 2023

Related U.S. Application Data
(65) Continuation of application No. 16/565,403, filed on Sep. 9, 2019, now Pat. No. 11,464,648.

ABSTRACT
A multi-portal method for treating a subject's spine includes distracting adjacent vertebrae using a distraction instrument positioned at a first entrance along the subject to enlarge an intervertebral space between the adjacent vertebrae. An interbody fusion implant can be delivered into the enlarged intervertebral space. The interbody fusion implant can be positioned directly between vertebral bodies of the adjacent vertebrae while endoscopically viewing the interbody fusion implant using an endoscopic instrument. The patient's spine can be visualized using endoscopic techniques to view, for example, the spine, tissue, instruments, and implants before, during, and after implantation, or the like. The visualization can help a physician throughout the surgical procedure to improve patient outcome.
(Continued)

Primary Examiner — Eduardo C Robert
Assistant Examiner — Steven J Cotroneo
(74) *Attorney, Agent, or Firm* — Perkins Coie LLP

Int. Cl. **A61B 1700** (2006.01)
A61B 1702 (2006.01)
(Continued)

U.S. Cl. **A61B 170224** (2013.01); **A61B 170218** (2013.01); **A61B 20170225** (2013.01)
(Continued)

Field of Classification Search
CPC **A61B 173421; A61B 117074**
See application file for complete search history.

References Cited
U.S. PATENT DOCUMENTS
3,486,505 A 12/1969 Morrison
5,171,279 A 12/1992 Mathews
5,209,021 A 2/1993 Pokuandi
5,437,017 A 8/1995 Lober et al.
5,699,232 A 11/1997 Tien
5,762,629 A 6/1998 Kambara
5,782,832 A 7/1998 Larson et al.
(Continued)

OTHER PUBLICATIONS
Choi, Chang Myung et al. "How I do it? Bipodal endoscopic spinal surgery (BESS) for treatment of lumbar spinal stenosis." *Acta Neurochir* (2016) 158:459-463, published Jan. 18, 2016.
(Continued)


Prior Publication Data
US 2023/029751 A1 Oct. 19, 2023

Related U.S. Application Data
(63) Continuation of application No. 16/687,520, filed on Nov. 18, 2019, now Pat. No. 11,678,906, which is a (Continued)

Int. Cl. **A61B 1734** (2006.01)
A61B 500 (2006.01)
(Continued)




U.S. Cl. **A61B 173423** (2013.01); **A61B 524** (2013.01); **A61B 50893** (2013.01);
(Continued)

Field of Classification Search
CPC **A61F 2444-447; A61F 24611; A61B 17025; A61B 20170256-0262; A61B 20173445; A61B 20173447**
See application file for complete search history.



U.S. Patent

Apr. 9, 2024 Sheet 5 of 13 US 11,950,770 B1



ABSTRACT
A multi-portal method for treating a subject's spine includes distracting adjacent vertebrae using a distraction instrument positioned at a first entrance along the subject to enlarge an intervertebral space between the adjacent vertebrae. An interbody fusion implant can be delivered into the enlarged intervertebral space. The interbody fusion implant can be positioned directly between vertebral bodies of the adjacent vertebrae while endoscopically viewing the interbody fusion implant using an endoscopic instrument. The patient's spine can be visualized using endoscopic techniques to view, for example, the spine, tissue, instruments, and implants before, during, and after implantation, or the like. The visualization can help a physician throughout the surgical procedure to improve patient outcome.
(Continued)

Primary Examiner — Eduardo C Robert
Assistant Examiner — Steven J Cotroneo
(74) *Attorney, Agent, or Firm* — Perkins Coie LLP

Int. Cl. **A61B 1700** (2006.01)
A61B 1702 (2006.01)
(Continued)

U.S. Cl. **A61B 170224** (2013.01); **A61B 170218** (2013.01); **A61B 20170225** (2013.01)
(Continued)

Field of Classification Search
CPC **A61B 173421; A61B 117074**
See application file for complete search history.

References Cited
U.S. PATENT DOCUMENTS
5,886,182 A 12/1998 Wolton
5,888,223 A 3/1999 Bupp
6,126,880 A 10/2000 Bupp
6,182,170 A 12/2000 Fisher et al.
6,348,038 B1 2/2002 Melker et al.
6,490,756 B1 6/2002 Brott (Continued)

FOREIGN PATENT DOCUMENTS
CN 105636555 A 6/2016 (Continued)
EP 3016617 A2 5/2016 (Continued)

OTHER PUBLICATIONS
Choi, Chang Myung et al. "How I do it? Bipodal endoscopic spinal surgery (BESS) for treatment of lumbar spinal stenosis." *Acta Neurochir* (2016) 158:459-463, published Jan. 18, 2016.
(Continued)

Prior Publication Data
US 2021/0068975 A1 Mar. 11, 2021






Related U.S. Application Data
(65) Continuation-in-part of application No. 16/565,403, filed on Sep. 9, 2019, now Pat. No. 11,464,648.

ABSTRACT
A multi-portal method for treating a subject's spine includes distracting adjacent vertebrae using a distraction instrument positioned at a first entrance along the subject to enlarge an intervertebral space between the adjacent vertebrae. An interbody fusion implant can be delivered into the enlarged intervertebral space. The interbody fusion implant can be positioned directly between vertebral bodies of the adjacent vertebrae while endoscopically viewing the interbody fusion implant using an endoscopic instrument. The patient's spine can be visualized using endoscopic techniques to view, for example, the spine, tissue, instruments, and implants before, during, and after implantation, or the like. The visualization can help a physician throughout the surgical procedure to improve patient outcome.
(Continued)


Primary Examiner — Andrew Yang
(74) *Attorney, Agent, or Firm* — Perkins Coie LLP

Int. Cl. **A61B 2444** (2006.01)
A61B 1702 (2006.01)
(Continued)

U.S. Cl. **A61F 24455** (2013.01); **A61B 10065** (2013.01); **A61B 17025** (2013.01); **A61B 20170256** (2013.01); **A61B 2017681** (2013.01)



The dualPortal Training Institute



The dualPortal Value

- Vast amounts of training and clinical reference resources including access to the biggest network of KOLs and the latest procedural innovations.
 - Annual Endo Symposium & dualPortal US TOUR
 - Regional labs – 60 labs in 2024
 - DTI – dualPortal Training Institute in Irvine, CA
 - Advance training trips to S. Korea
 - KOL-surgery observation/mentorship programs in US and Korea
 - Latest technique updates from the experts in Korea and US
- The most advanced technologies available in the US with easy to use products that gets updated for free as they become available in the future
- The most dependable and continuous support system both from the HQ and local sales/tech representatives.



Cost effectiveness of dualPortal






1. Overall, hospital savings via less LOS, less complications, higher patient satisfaction (publication below)
2. More patient flow via marketing of endo solutions
3. Flexible ways to maximize discount with Amplify
 - Rebate systems tied to implant/dualX cage usage
 - Volume discount
 - Endo fusion related discounts
 - Discounts on disposables tied to one time capital purchase

<https://thejns.org/spine/view/journals/j-neurosurg-spine/40/1/article-p77.xml>

1. Endoscopic procedures were more costly for hospitals, with disposable supplies accounting for 31.7% of the total cost of cases, compared to 10.1% of the cost for open procedures.
2. A multivariate linear regression independently associated endoscopic procedures with higher total costs.
3. Patients who underwent the open procedures had longer hospital stays, at 1.4 days on average, compared to 0.7 day for patients who underwent the endoscopic procedure.
4. Nearly 8% of patients who underwent the open procedure reported perioperative complications, compared to 3.1% of the endoscopic patients.

•

The Cage Dilemma




ALIF



TLIF

dualX: Designed to Solve Key Challenges of Expandable Interbody Devices



High Subsidence / Collapse Risk with Uni-directional Expandable Cages

The incidence of cage subsidence was higher in the expandable group (19.7% vs 5.4%, $p = 0.0017$). Within the expandable group, the unilateral facetectomy-only subgroup had a 5.6 times higher subsidence rate than the PCO subgroup (26.8% vs 4.8%, $p = 0.04$). Four expandable cages collapsed over time.

dualX Broad Footprint Mitigates the Risk of Subsidence

Long-term radiographic outcomes of expandable versus static cages in transforaminal lumbar interbody fusion

Chih-Chang Chang^{1, 2, 3, 4}, Dean Chou¹, Brenton Pennicooke¹, Joshua Rivera⁵, Lee A Tan¹, Sigurd Berven⁶, Praveen V Mummaneni¹

Affiliations + expand

PMID: 33186902 DOI: 10.3171/2020.6.SPINE191378

Abstract

Objective: Potential advantages of using expandable versus static cages during transforaminal lumbar interbody fusion (TLIF) are not fully established. The authors aimed to compare the long-term radiographic outcomes of expandable versus static TLIF cages.

Methods: A retrospective review of 1- and 2-level TLIFs over a 10-year period with expandable and static cages was performed at the University of California, San Francisco. Patients with posterior column osteotomy (PCO) were subdivided. Fusion assessment, cage subsidence, anterior and posterior disc height, foraminal dimensions, pelvic incidence (PI), segmental lordosis (SL), lumbar lordosis (LL), pelvic incidence-lumbar lordosis mismatch (PI-LL), pelvic tilt (PT), sacral slope (SS), and sagittal vertical axis (SVA) were assessed.

Results: A consecutive series of 178 patients (with a total of 210 levels) who underwent TLIF using either static (148 levels) or expandable cages (62 levels) was reviewed. The mean patient age was 60.3 ± 11.5 years and 62.8 ± 14.1 years for the static and expandable cage groups, respectively. The mean follow-up was 42.9 ± 29.4 months for the static cage group and 27.6 ± 14.1 months for the expandable cage group. Within the 1-level TLIF group, the SL and PI-LL improved with statistical significance regardless of whether PCO was performed; however, the static group with PCOs also had statistically significant improvement in LL and SVA. The expandable cage with PCO subgroup had significant improvement in SL only. All of the foraminal parameters improved with statistical significance, regardless of the type of cages used, however, the expandable cage group had greater improvement in disc height restoration. The incidence of cage subsidence was higher in the expandable group (19.7% vs 5.4%, $p = 0.0017$). Within the expandable group, the unilateral facetectomy-only subgroup had a 5.6 times higher subsidence rate than the PCO subgroup (26.8% vs 4.8%, $p = 0.04$). Four expandable cages collapsed over time.

Conclusions: Expandable TLIF cages may initially restore disc height better than static cages, but they also have higher rates of subsidence. Unilateral facetectomy alone may result in more subsidence with expandable cages than using bilateral PCO, potentially because of insufficient facet release. Although expandable cages may have more power to induce lordosis and restore disc height than static cages, subsidence and endplate violation may negate any significant gains compared to static cages.

Post Operative Collapse/Loss of Height Uni-directional Expandable Cages

JNS JOURNAL OF
NEUROSURGERY
OFFICIAL JOURNALS OF THE AANS SINCE 1944

Journal of Neurosurgical Spine, 2020 Nov 13: 1-10





FIG. 3. Lateral radiographs 1 month after TLIF with an expandable cage (A), 12 months after TLIF with an expandable cage demonstrating cage collapse/loss of height of the cage (B), 40 months after TLIF demonstrating pseudarthrosis and screw fracture (*arrow*; C), and after revision surgery with anterior lumbar interbody fusion as well as posterior screw-rod replacement, demonstrating solid fusion (D).

dualX TLIF cage: A Revolution in Expandable Interbody Devices

Minimize Subsidence –
Wide Horizontal Expansion
Largest Footprint

trueLordosis™
8°, 12°, 15°, 18°*



**Long Term Durability,
Stability**
– Two Independent Locking
Mechanisms

**Highest
Post-Expansion Graft
Volume Delivery**

AMPLIFY™

SURGICAL





TRANSFORMING THE
→
ORDINARY



dualX – The Largest Footprint Expandable Cage



Safe and Secure

- Minimize subsidence due to wide footprint
- Only implant that provides wide horizontal expansion followed by powerful vertical expansion
- Allows for completely endoscopic placement





"...like you're performing an ALIF through the back!"

ALIF-like TLIF



trueLordosis™

- The Implant will always provide the prescribed Lordosis
- Available in 8°, 12°, 15°, 18°*



dualX – Long Term Durability and Inherent Stability

Ensures Durability and Stability with Two Independent Locking Mechanisms


dual Locking

Expansion Locking Mechanism



Secondary Screw Lockout

- Innovative dual locking design
 - Maintains the integrity of the implant until the patient is fused
- Final locking screw
 - Ensures implant stays expanded in width and height
- The only one of two “non-screw based” Expansion Mechanism



Post-expansion, Surgeon Preferred Bone Grafting

Maximize Bone Graft Delivery

- Integrated Post Packing Through Delivery Handle
- Large Internal Atrium Retains Extensive Bone Graft Volume
- Unique “Open Structure” Enables Bone Graft to Flow Beyond Cage and Fill Entire Disc Space





dualX[®]Slim


AMPLIFY[®]
SURGICAL

dualX[®]Slim


dualPortal Fusion / dualLIF



dualPortal Fusion / dualLIF


































dualLIF®



dualPortal® + dualX®
dual[LIF]

Market Leading, Differentiated Benefits

		Accelus	Globus	Medtronic	Nuvasive
Bi-Directional Expansion					
Large Footprint					
Largest Footprint Size (WxL) (vs. height expanding devices)	21x22/27mm	14x29	12x30	10x32	11x36
Significant Volume for Internal Bone Graft Filling					
Dual Locking Safety					
Solution to Minimize Psoas Retraction (LLIFs)					
All Titanium Solution and Adaptable for 3D Printing					

dualX LLIF with Small Portal Access

dualX[®]

Dual Expanding Interbody Fusion System

LLIF

LATERAL LUMBAR INTERBODY FUSION

dualX[®] transforms the fusion environment from insertion to spinal restoration by delivering a powerful **dual-expanding** implant through a minimally invasive approach.

SMALL PORTAL ACCESS SYSTEM

Designed to minimize neural retraction while still achieving maximum implant geometry (22mm width, expanded).



dualX[®] LLIF


- Wide, horizontal expansion engineered to reduce subsidence risk.
- Powerful vertical expansion restores disc height for decompression.
- Highest amount of post-expansion, surgeon-preferred bone grafting.
- Innovative dual-locking mechanism designed to maintain integrity of implant.



18MM OPENING

The procedure is performed completely through an **18mm** opening, reducing potential psoas retraction up to 50%.

LOCKING SCREW



Available at 0°, 7°, and 12° lordotic angles for sagittal alignment restoration.

AMPLIFY[®]
SURGICAL



dualX – A New LLIF Leader

Today...



<40% of Surgeons Practice LLIF

- Want to avoid post-op complications
 - 30% incidence
- Overstretching psoas
 - 30mm+ retraction required
- No other devices solve issue
- Nuvasive dominates the market

dualX LLIF Solution


Eliminates excessive psoas retraction

- **Safer:** Minimize thigh related complications
- **Faster:** Simpler set-up, less intra-op X-rays
- **Economic:** Low-cost retractor tube



dualX ALIF / ATP: FDA Cleared


Anterior Dual-Expanding Interbody



Designed to Minimize
Retraction of Major Vessels (L3/4, L4/5)

FDA
Cleared

ATP Dual-Expanding Interbody




- ATP-specific Lordosis
- Only Cage Specifically Shaped for ATP
 - Safer access through dangerous anatomy
 - Application for endoscopic access
- Obviates Need to Conduct Orthogonal Sweep
 - Less trauma to the Psoas



Expanding the Opportunity for MISS Surgery


dualX[®]
Dual Expanding Interbody Fusion System




Instrument Simplicity & Safety

- All steps performed safely through a single inserter
 - Insertion
 - Lateral expansion
 - Vertical expansion
 - Graft filling
 - Screw lock out


Common, Familiar
Inserter for Every
Surgical Approach



One Instrument Set for All Five dualX Systems



Simply Swap Top Tray




Surgeon Experience the Benefits of dualX

"I have done over 60 levels with the entire range of Amplify cages.... I have not had to revise any and all that are beyond 6 months post op are considered fused without any signs of non-union.

- Surgeon from AZ

I've used about 100 cages so far with lots of follow-ups past 12-18+ months. The footprint is impressive. Much more stable and haven't seen any subsidence. The larger footprint helps, and more stability helps with post-op pain. It also provides impressive power to correct deformity.

- Surgeon from IA



Building Clinical Momentum

- Operating in over 25 states



March 29 | UCI Health | 7:00 AM

AMPLIFY[®] SURGICAL 4th Annual Endoscopic Spine Symposium
featuring dualPortal & dualX technologies

Register: <https://amplifysurgical.com/4th-annual-symposium/>

Symposium Faculty



Symposium Chair
Dr. Don Y. Park
Orange, CA



Dr. Cheol Woong Park
Daejeon, S. Korea



Dr. Dong Hwa Heo
Seoul, S. Korea



Dr. Charla Fischer
New York, NY



Dr. Samuel Cho
New York, NY



Dr. Jon J.W. Yoon
Philadelphia, PA



Dr. Nam Lee
Busan, S. Korea



Dr. Man Kyu Park
Busan, S. Korea



Dr. Andrew Chung
Glendale, AZ



Dr. Ki Eun Chang
San Diego, CA



Dr. Gregory Basil
Miami, FL



Dr. Sohaib Hashmi
Orange, CA



Dr. Young San Ko
Daegu, S. Korea

Questions: marketing@amplifysurgical.com

AMPLIFY


SURGICAL

4th Annual Amplify Surgical Endoscopic Symposium


Featuring dualPortal & dualX technologies

SAT. MARCH 29, 2025

Register Here:
amplifysurgical.com/4th-annual-symposium/



dualPortal



*Hands on
Cadaver Workshop
A multi-station and
multi-session lab course*



Dr. Don Park, 500th dualPortal Procedure

- “Over the past 100 cases, I have experienced how dualPortal spinal endoscopy can push the boundaries of spine surgery to a whole different level. My practice has completely transformed due to the enhanced visualization and the shift to the outpatient ASC setting. From endoscopic discectomies to fusions, my patients have benefited so much from this enabling technique. dualPortal spinal endoscopy is the next evolution of minimally invasive spine surgery.”*

CONGRATULATIONS

Dr. Don Y. Park



on completing **500**
dualPortal[®] CASES!



Dr. Don Y. Park
 UCI Health
 Orange, CA

Endoscopic Approach

500



dualPortal® Clinical Paper

European Spine Journal
<https://doi.org/10.1007/s00586-023-07701-9>

REVIEW ARTICLE

Clinical outcomes and complications after biportal endoscopic spine surgery: a comprehensive systematic review and meta-analysis of 3673 cases

Don Y. Park¹ · Alexander Uphill-Brown¹ · Nora Curtin¹ · Christopher D. Hamad¹ · Akash Shah¹ · Brian Kwon² · Yong H. Kim³ · Dong Hwa Heo⁴ · Cheol Woong Park⁵ · William L. Sheppard¹

Received: 7 December 2022 / Revised: 7 December 2022 / Accepted: 4 April 2023
 This is a U.S. Government work and not under copyright protection in the US; foreign copyright protection may apply 2023

Table 2 Operative Detail

Biportal endoscopic discectomy	Biportal endoscopic decompression	Biportal endoscopic TLIF	Surgery performed at lumbar levels				
			L1-L2	L2-L3	L3-L4	L4-L5	L5-S1
1056	2402	261	2	101	351	1255	338

Table depicting the number of each surgery type performed and the lumbar levels involved in the surgeries

Our study examining **3673 cases** demonstrated that **biportal spinal endoscopy is safe and effective**. Complications are **comparable** or **less** than previously published rates with respect to open and microscopic techniques...**Significant improvement in VAS-Back, VAS-Leg, ODI, and Macnab Scores** were seen across the cohort.

AMPLIFY
SURGICAL

Clinical Study Indicates dualPortal® Spinal Endoscopy as Successful Surgical Technique Demonstrating Favorable Patient Outcomes and Low Complication Profile

Recent clinical study indicates positive data on Amplify Surgical's dualPortal® Endoscopic Approach.

IRVINE, CA (PRWEB) - April 28, 2023

Amplify Surgical, Inc., a medical device company focused on innovative minimally-invasive surgery for the lumbar spine, announced positive data from the first comprehensive systematic review and meta-analysis on its dualPortal® Spinal Endoscopy, also known clinically as 'biportal spinal endoscopy'. The study, "Clinical Outcomes and Complications after Biportal Endoscopic Spine Surgery: A Comprehensive Systematic Review and Meta-Analysis of 3673 Cases," was recently published in the *European Spine Journal*, led by principal investigator, Don Y. Park, M.D., Vice Chair of Quality and Safety of the Department of Orthopaedic Surgery, David Geffen School of Medicine at UCLA in Los Angeles, California. Link: <https://pubmed.ncbi.nlm.nih.gov/37079079/> and <https://amplifysurgical.com/news/press-release-23/>.

AMPLIFY
SURGICAL



dualLIF® Clinical Paper

Neurospine 2023;20(1):110-118.
<https://doi.org/10.14245/ns.2346116.058>

Neurospine
 pISSN 2386-6383 eISSN 2386-6391

Original Article

The Use of Dual Direction Expandable Titanium Cage With Biportal Endoscopic Transforaminal Lumbar Interbody Fusion: A Technical Consideration With Preliminary Results

Corresponding Author
 Dong Hwa Heo
<https://orcid.org/0000-0003-1203-4550>

Department of Neurosurgery, Endoscopic Spine Surgery Center, Champodonamu Hospital, 32 Baumae-ro 35-gil, Seocho-gu, Seoul 06744, Korea
 Email: spineurosurgeon@gmail.com

Received: January 25, 2023
 Revised: February 27, 2023
 Accepted: February 27, 2023

Don Young Park¹, Dong Hwa Heo²

¹Department of Orthopaedic Surgery, David Geffen School of Medicine at UCLA, Los Angeles, CA, USA
²Department of Neurosurgery, Spine Center, Champodonamu Spine Hospital, Seoul, Korea

Table 2. Clinical results

Variable	Preoperative	Postoperative		
		6 Weeks	3 Months	6 Months
VAS back*	6.9 ± 1.19	2.1 ± 1.85	1.3 ± 1.57	1.25 ± 0.63
VAS leg*	8.3 ± 1.16	0.55 ± 1.57	1.6 ± 1.65	1.0 ± 0.94
ODI*	55.2 ± 9.1	32.3 ± 17.3	29.1 ± 15.5	26.6 ± 7.5

Values are presented as mean ± standard deviation.
 VAS, visual analogue scale; ODI, Oswestry Disability Index.
 *p < 0.05.

"There was no significant subsidence or collapse of the expandable cages during the 6-month follow-up period. Lumbar lordosis and disc height were significantly increased after surgery. ODI and VAS scores were significantly improved at 6 months after surgery."

AMPLIFY
SURGICAL

Clinical Study Demonstrates Favorable Patient Outcomes with the dualPortal® Endoscopic TLIF with dualX®, the dualLIF® System

Recent clinical study indicates positive data on Amplify Surgical's dualLIF® System, dualPortal® Endoscopic TLIF with dualX®.

IRVINE, CA (PRWEB) – April 14, 2023

Amplify Surgical, Inc., a medical device company focused on innovative minimally-invasive surgery for the lumbar spine, today announced positive data from a retrospective study on its dualLIF® System, the synthesis of dualPortal® Spinal Endoscopy and dualX® Dual-Expanding Interbody Technology. The study, "The Use of Dual Direction Expandable Titanium Cage With Biportal Endoscopic Transforaminal Lumbar Interbody Fusion: A Technical Consideration With Preliminary Results," was recently published in *Neurospine*, led by clinical researchers, Don Y. Park, M.D., Vice Chair of Quality and Safety of the Department of Orthopaedic Surgery, David Geffen School of Medicine at UCLA in Los Angeles, California and Dong Hwa Heo, M.D., Ph.D., Director of the Department of Neurosurgery, Spine Center, Champodonamu Spine Hospital in Seoul, South Korea.

Link: <https://pubmed.ncbi.nlm.nih.gov/37016859/>

dualPortal®

Endoscopic Approach



Shallow
Learning Curve



Greater
Flexibility



Enhanced
Visualization




dualX Slim®


Dual Expanding
Interbody Fusion



#easyInsertion



#wideFootprint



#trueLordosis™
8° 12° 15°



#oneStepExpansion



Transforming the Ordinary